

Superior Craftsmanship



VETERINARY INSTRUMENTS CATALOGUE



Ask For it!

info@bergmed.com

Our Products Are Guaranteed And Meet The Relevant
Worldwide Standards And Our Products Are
Manufactured Entirely Under In Factory.

Company profile



BERGMED International is a leading manufacturer and supplier of all kinds of Veterinary, Surgical & Dental Instruments, T.C, and Orthodontic tools, Implants, Hollow-wear and wide range of Scissors, Pliers & Cutters to worldwide. All of our instruments are ISO 9001:2008, ISO 13485:2018 and CE mark certified in accordance with the Medical Devices Directive. We also provide OEM services and customized packaging to our valuable customers. We offer special discounted prices for our distributors.

Having our manufacturing unit in Sialkot (Pakistan), world's renowned city for the manufacturing of Veterinary, Surgical & Dental Instruments. We worked hard to build our reputation in this evolving industry and our satisfied customers are our strength. Our products are manufactured using the best quality stainless steel as AISI, ASTM and DIN standard's recommended which is imported from Germany, Japan and France. The key to our success is our quality assurance. We use the right material for the right products.

Health and Safety Policy, Environment Policy always adopted as per international standard.

Quality Policy:-

- Improved Quality
- Cost Reductions
- Working for Customer satisfaction
- Customized products for special needs customers
- Package designing and innovation
- New products or existing product enhancements

Bergmed International prices are very competitive and very flexible with customers provide maximum co-operation. We aim to provide excellent customer service and we promote a 'can do' culture. We aspire to being the easiest company in the market to deal with for any and all of your surgical instrument's needs.

Our customers are our most important asset and **BERGMED International** endeavors to provide a personalized customer service experience. We aim to provide the highest level of customer service to the International Healthcare market and are proud that we have built and continue to maintain solid relationships with our customers and distributors.

At all times we are professional, honest, ethical and reliable.

Thanks and Best Regards,

CONTENTS



01	Balling Guns	04
02	Bone Surgery	05-06
03	Castration Instruments	07-08
04	Ear Marking & Taging	09
05	Dehorner Instruments	10
06	Equine Dental Instruments	11
07	Farrier Tools	12-14
08	Grooming Equipment & Hoof Rasps	15
09	Hoof & Claw Instruments	16-20
10	Injectors & Extractors	21
11	Obstetrical Instruments	22
12	Plaster Instruments & Snake Catchers	23
13	Restraining Instruments	24
14	Ring Applicators & Teat Instruments	25
15	Vaginal Speculums	26
16	Snaffles. Dee. Ring Snaffle Bits	27
17	Spurs and Stirrups	29
18	Surgical Forceps & Towel Clamps	30
19	Dressing and Tissue Forceps, Hollowares, Scalpel Handles	31
20	Needle Holders & Scissors	32
21	Retractors	33
22	Equine Surgical & Miscellaneous Instruments	34
23	Root Elevators	35
24	Scalers, Probes, Dental Extraction Forceps & Mouth Retractors for Cats & Dogs	37

Balling Guns

BRI-1-1001 Mini Balling Gun

Balling gun with spring clips 1.25" Metal head, brass : Length 14.75"



BRI-1-1002 Metal Head Balling Gun

Balling gun with 3" Metal head, brass: Length 17".



BRI-1-1003 Plastic Head Balling Gun

Balling gun with 2.5" Plastic head, brass: Length 17".



BRI-1-1004 Balling Gun

Balling gun with 5" Metal head, brass: Length 22".



BRI-1-1005 Rubber Head Balling Gun

Balling gun with spring clips 4" Rubber head, brass: Length 23.5".



BRI-1-1006 Pistol Balling Gun / Multi Bolus Balling Gun


Nickel-plated metal balling gun which holds up to six large boluses at one time. Uses a spring-action end-clip to secure boluses in barrel: Length 18.25".



BRI-1-1007 British Capsule Forceps

Excellent for capsuling sheep or horses. Stainless steel: Length 12".






BRI-1-1008 Bone Cutting Forceps (Stille-Liston)

Also available in curved blade. Length 10.5".




BRI-1-1009 Bone Cutting Forceps (Stille-Liston)

Angled blades. Length 10.5"



BRI-1-1010 Surgical Wire Cutters

Forged from stainless steel and reinforced with tungsten carbide tip or specially



BRI-1-1011 Surgical Wire Cutters

Forged from stainless steel and reinforced with tungsten carbide tip or specially tempered tool steel. Size 9.25"




BRI-1-1012 Bone Rongeur

(Leksell-Stille) Length 9" s/s.




BRI-1-1013 Surgical Pin Cutters

Cutting tools used to cut the wires and pins implanted during the orthopedic surgeries to fix the fractures. They are forged from stainless steel and reinforced with tungsten carbide or specially tempered tool steel. The handles are so firm that they are able to provide more power for cutting and absorb the jerk of wire cutting. Length 18.5".



BRI-1-1014 Bone Rongeur Forceps (Stille)

Angled on side. 6x12mm jaws. Length 9" s/s.



BRI-1-1015 Bone Holding Forceps (Verbrugge)

Verbrugge bone forceps for holding / positioning re-positioning bone. Available sizes 10.5" s/s.



BRI-1-1016 Bone Holding Forceps

Length 11" s/s



BRI-1-1017 Bone Holding Forceps (Lambotte)

Lambotte bone holding forceps with ratchet. Available Sizes 8.5", 10", 13" s/s.

Bone Surgery



BRI-1-1018 Bone Holding Forceps
Length 13.5" s/s.



BRI-1-1019 Bone Holding Clamp Lowman
Lowman bone clamp for holding, positioning, re-positioning bone. 1 x 2 Prong Jaws. Length 8" s/s.

Branding Instruments



BRI-1-1020 Revolving Head Tattoo Pliers
Made of special aluminum alloy. Head holds four dies on each side, for a total of eight digits. Size 3/8" revolving head tattoo outfit.



BRI-1-1021 Standard Tattoo Pliers
Designed in the same general style as our regular cattle tattoos but reduced in size to make work more convenient with small animals. Head holds four dies, for a total of four digits. Size 3/8".



BRI-1-1022 Small Tattoo Pliers
Light weight pattern, designed for small animals. Made of aluminum alloy. Head holds four dies, for a total of four digits. Size 0.25".



BRI-1-1023 Tattoo Letter Set 0-9
Tattoo digits for Revolving head tattoo pliers and Standard tattoo pliers. Size 3/8". Tattoo digits for Small tattoo pliers. Size 0.25".



BRI-1-1024 Aluminum Paint Branding Irons
Made of aluminum castings, light in weight. Set consists of nine irons in 2.5" sizes. 0.5" depth, 1.5" width. Also available in 4" sizes. NOTE: A set 0-9 does not include the figure 9. The 6 iron can be reversed to make the figure 9.



BRI-1-1025 Aluminum Tattoo Pliers
The Aluminum Tattoo Pliers has a heavy duty spring for quicker release.



BRI-1-1026 Tattoo Pliers
Made of carbon steel.



BRI-1-1027 Tattoo Pliers
Two rows of tattoo digits can be inserted at a time. Made of carbon steel.

Castration Instruments



BRI-1-1028 Emasculatome 9" Castrator- Extra Small Size

These are precision-made instruments for the humane and bloodless castration of live stock by severing the testicular cord without injury to the bag (scrotum). Also Available with plastic and wooden grips. For lambs. Length 9".



BRI-1-1029 Emasculatome 12" Castrator- Small Size

Also available with plastic grips. For lambs and young calves or for docking lambs tail. Length 12".



BRI-1-1030 Emasculatome 14" Castrator- Medium Size

Also available with plastic grips. For calves and young bulls. Length 14".



BRI-1-1031 Emasculatome 16" Castrator- Large Size

Also available with plastic grips. For mature bulls. Length 16"



BRI-1-1032 Emasculatome 19" Castrator- X Large Size

Also available with plastic grips. For bulls, horses, buffalo, and camels. Length 19".



BRI-1-1033 Large/Deluxe Castrating Band Applicator Aluminum

Heavy-duty aluminum ring applicator with notched expanding fingers to hold the ring. Useful on lambs, calves, goats and other animals



BRI-1-1034 Castrating Band Applicator-Stainless Steel

For applying latex rings/bands used in tail docking and castration. The tapered nose prevents dirt build-up. Roller insures smooth operation: Length 8.75".



BRI-1-1035 Castrating Band Applicator- Nickel Plated

For applying latex rings/bands used in tail docking and castration. The tapered nose prevents dirt build-up. Roller insures smooth operation: Length 8.75".



BRI-1-1036 Blue Economy Castrating Band Applicator

Designed for quick and easy one-person operation. Lightweight blue plastic grips, textured for sure-grip, comfortable to use. Useful on lambs, calves, goats and other animals.



BRI-1-1037 Crutch Knee Wrench - For Large Emasculatomes

The use of the crutch knee attachment makes possible a one-person operation with the larger emasculatome. Set consists of polished crutch attachment and spanner wrench. Fits 14", 16" and 19" emasculatomes only.

Castration Instruments



BRI-1-1038 Ear Docking

For dog's ear docking. Aluminum made. Straight, 12cm. Also available in curved.



BRI-1-1039 Frank's Emasculator

Length 7-1/2", 9", 10-1/2" s/s.



BRI-1-1040 Haussman's Emasculator

Double crush prevents primary and secondary hemorrhage because of the multiple crushing of the spermatic cord and artery. Sizes available 7-1/2", 9", 10", 11", 12" s/s



BRI-1-1041 White's Emasculator

Length 7-1/2", 9", 10", 11", 12" s/s.



BRI-1-1042 Haussman's Emasculator / Real Product

Curved, sizes available from 7-1/2", 9", 10, 11", 12" s/s.



BRI-1-1043 Reimer's Emasculator

A very popular emasculator has a double crushing part. The cutting lever is worked separately. This allows extra time for crushing action and eliminates after bleeding: Length 12" s/s.



BRI-1-1044 Serra's Emasculator Curved

A concave cutting plate draws the cord concentrically in to the jaws of the castrator. Reduces tissue slippage and post-op hemorrhage. Length 14" s/s



BRI-1-1045 Serra's Emasculator Curved With Ratchet

It has curved non slip grips and positive crush for one hand use with ratchet: Length 14" s/s.



BRI-1-1046 Serra's Emasculator With Ratchet

This newly designed instrument draws the cord concentrically into the jaws, thus preventing slipping of the cord. ratchet holds emasculator locked: Length 13-1/2" s/s.



BRI-1-1047 Sand's Emasculator

Draws the cord concentrically in to the jaws, thus preventing slipping of the cord. Hemorrhage is reduced to a minimum ratchet holds forceps locked: Length 12" s/s



BRI-1-1049 Gas Tail Docker

For quick, painless and economical tail docking of piglets. Stainless steel.



BRI-1-1048 Castration Hook for Male cats



BRI-1-1050 Newberry Castrating Knife

A stainless steel castrating knife will split the scrotum for easier removal of testicles with 2" blades. Overall Length 10.5".

Ear Marking & Taging



BRI-1-1051 Flower-Cut Ear Notcher

The stainless steel construction with sharp, clean cut and slide spring action handle make these notchers ideal for identification purposes. Cuts notch 1" long. Length 8".



BRI-1-1052 W-Cut Ear Notcher

A high-quality, stainless steel ear notcher Cuts a "W" notch 0.70" long. Length 6".



BRI-1-1053 U-Cut Ear Notcher

A high-quality, stainless steel ear notcher Cuts a "U" notch 0.70" long. Length 6".



BRI-1-1054 Small Ear Notcher With V-Cut

A high-quality, stainless steel ear notcher Cuts a "V" notch 8mm/0.32" base x 16mm/0.63" long. Length 6".



BRI-1-1055 Medium Ear Notcher With V Shape Die

A high-quality, stainless steel ear notcher Cuts a "V" notch 12mm/0.47" base x 18mm/0.71" long. Length 6".



BRI-1-1056 Large Ear Notcher With V Shape Die

carefully ground after being hardened and tempered. Dies are not interchangeable. 20.85mm/0.82" cut. Length 9.37".



BRI-1-1057 Ear Tag Applicator



BRI-1-1058 Plastic Tag Applicator

For applying plastic tags to the ears. Length 9".

Dehorers Instruments



BRI-1-1059 Dehorner Wooden Handle
Wooden handle dehorner: Length 13" c/s.



BRI-1-1060 Dehorner Metal Handle
Straight metal handle dehorner with rubber grips: Length 13" c/s, c/p.



BRI-1-1061 Dehorner Metal Handle
Straight metal handle dehorner with rubber grips: Length 17" c/s, c/p.



BRI-1-1062 Dehorner Metal Handle
Straight metal handle dehorner with rubber grips. Heat-treated steel cutting blades: Length 13" c/s, c/p.

Equine Dental Instruments



BRI-1-1063 Tooth Cutter
Hardened stainless steel simple joint lower molar cutters with permanent handles. Length 22" s/s.




BRI-1-1064 Tooth Cutter
Molar cutter with simple joint. Available in 3 different models. 1) Closed Edges. / 2) Medium Gap. / 3) Wide Gap. Length 22" s/s.



BRI-1-1065 Molar Forceps
Length 18", 22" s/s.



**BRI-1-1066 Equine Molar Forceps-
Gunther Tooth Forceps.** Length 22" s/s.



BRI-1-1067 Molar Forceps
Serrated jaw creates good purchase on cheek teeth and premolar caps. A hugely versatile pair of forceps for every day use. Length 22" s/s.



BRI-1-1068 Molar Spreader
Used for loosening cheek teeth prior to extraction: Length 19", 22" s/s.



BRI-1-1069 Molar Incisor Extractor
These are designed to remove hard to get at rear molars but can also be used to remove incisors. Length 22" s/s.



BRI-1-1070 Equine Compound Molar Cutter
Stainless steel compound molar cutters with permanent handles. Powerful compound action. Available in 4 different models. 1) Fully Open Jaw. / 2) Closed Jaw. 3) Part Open Jaw. / 4) Part Open Jaw Angled. Length 24", 25", 26" s/s.

Equine Dental Instruments



BRI-1-1071 Wolf Tooth Forceps

These extracting forceps are ideal for use in wolf teeth extractions and incisor cap extractions. Length 9" s/s.



BRI-1-1073 Cap Extractor

Lower (Renold's) Length 17" s/s



BRI-1-1075 Compound Rongeurs

Stainless steel compound rongeur with offset head works well for wolf teeth, blind wolf teeth, incisor fragments and most tooth and cap slivers. Length 9" s/s.



BRI-1-1077 3 Root Box Joint Forceps On Side

With knurled handles. For molars with very little crown in aged horses and deciduous premolar caps: Length 19" s/s



BRI-1-1079 Dental Float Moveable

Float comes with the moveable head with one wooden handle. With two screws in back to hold blade. Length 22".



BRI-1-1081 Dental Float

Combination float comes with both a straight and angled head. Length 10", with one wooden handle. Length 10.5". With two screws in back to hold blade. One rasp/rasp and one rasp/file in set. Length 20.5"



BRI-1-1072 Wolf Tooth Forceps / Extracting Forceps

Knurled handles provide optimum grip even when wet. Length 10" s/s.



BRI-1-1074 Cap Extractor

Upper (Renold's) Length 17" s/s.



BRI-1-1076 Fulcrum Serrated Jaw Box Joint Forceps

For extraction of first and second premolars or the number six and seven cheek teeth. This unit uses the hinge portion of the forceps for a fulcrum.



BRI-1-1078 3 Root Box Joint Forceps Off Side

With knurled handles. For molars with very little crown in aged horses and deciduous premolar caps. Length 19" s/s.



BRI-1-1080 Pig Tooth Nipper

Hardened stainless steel nippers with no gap in the jaws, construction with spring handle action. Also for cutting tails. Length 5" s/s.



BRI-1-1082 Equine Speculum By Schoupe

Reduces the danger of slipping with this instrument since it is held in place by the equine cheek area. I



BRI-1-1083 Bayer's Mouth Wedge

Very handy when examining the oral cavity or doing dental work in equine, s/s.



BRI-1-1084 Hoof Nipper T.C Tip Modified

A modified instrument which is forged from stainless steel and reinforced with tungsten carbide tip. Length 15". Also available in high carbon steel.



BRI-1-1085 Hoof Nipper

Particularly handy for farriers trimming larger feet or cob type horses with thicker walls. Forged from stainless steel. Length 15". Also available in high carbon steel.



BRI-1-1086 Hoof Nipper

Ideal for trimming away excess hoof wall or claws. Forged from stainless steel. Length 14". Also available in high carbon steel.



BRI-1-1087 Hoof Nipper

A very nice precision instrument, forged from stainless steel. Length 12".



BRI-1-1088 Hoof Nipper

Length 14". Also available in 12", 15" c/s.



BRI-1-1089 Hoof Nipper

Length 12". Also available in 14", 15" c/s.



BRI-1-1090 Hoof Nipper Curved Blade T.C Tip Modified

A modified instrument which is forged from stainless steel and reinforced with tungsten carbide tip. Length 15". Also available in 12", 14" c/s.



BRI-1-1091 Hoof Nipper Curved Blade

A modified instrument which is forged from stainless steel and has a curved blade. Length 15". Also available in 12", 14" c/s.



BRI-1-1092 Heavy Duty Hoof Nipper

Fine quality steel hoof nipper with cushion grip handles. Length 15". Also available in 12", 14" c/s.



BRI-1-1093 Hoof Nipper

Designed exclusively to trim away excess hoof wall. Length 15". Also available in 12", 14" c/s.



BRI-1-1094 Hoof Nipper

Hoof nipper is used to trim away excess hoof wall. Vinyl handles give a comfortable feel in the hand. Length 15" c/s.



BRI-1-1095 Hoof Nipper

Hoof Nipper Length 12" c/s.



BRI-1-1096 Shoe Puller And Spreader With T.C Tip Modified

A combination shoe puller and spreader. Heat treated for strength. It is a modified instrument with a Tungsten Carbide tip. Length 13".



BRI-1-1097 Shoe Puller And Spreader

A combination shoe puller and spreader. Has exclusive pattern with sharp teeth on outside edges. Length 13" c/s.



BRI-1-1098 Shoe Puller And Spreader

Has exclusive pattern with sharp teeth on outside edges. For ease in spreading most shoe sizes and styles. Length 13" c/s.



BRI-1-1099 Shoe Pullers (Pull-Offs)

As sharp as a nipper so you can easily use it as a nail cutter. Length 15".

Farrier Tools



BRI-1-1100 Gooseneck Clincher
Length 12" c/s. Also available in stainless steel.



BRI-1-1102 Farrier Clincher
For use on high and low nails. Used to roll over nails with curved jaw. Length 13" c/s. Also available in stainless steel.



BRI-1-1104 Farrier Clincher
For use with larger nails used in saddle horse hoofs to roll over and clinch nails. Oil quoting makes it rust resistant. Length 12" c/s. Also available in stainless steel.



BRI-1-1106 Farrier Clincher
For use with large nails used in saddles horse hoofs to roll over and clinch nails. Oil quoting makes it rust resistant. Length 12" c/s. Also available in stainless steel.



BRI-1-1108 Cutting Nipper
For cutting and pulling with cushion grip handles. Length 12" c/s. Also available in stainless steel.



BRI-1-1110 Crease Nail Puller powder
Coated Removes nails from crease of shoes that has been on horse for some time, 12" carbon steel powder coated. Also available in stainless steel.



BRI-1-1112 Clinch Cutter
Clinch cutter on one end and pritchel on opposite end s/s.



BRI-1-1101 Low Nail Curved Jaw Clincher
The curved jaw clincher features a proven design for excellent performance and reliability. Length 13" c/s. Also available in stainless steel.



BRI-1-1103 Gooseneck Clincher
Powder coating makes it rust resistant. Length 12" c/s. Also available in stainless steel.



BRI-1-1105 Farrier Clincher
To roll over nails with curved jaw. Red vinyl coated handle gives a comfortable grip. Length 13" c/s. Also available in stainless steel.



BRI-1-1107 Cutting Nipper
For cutting and pulling with cushion grip handles. Oil tempered and drop forged. Length 12" c/s. Also available in stainless steel.



BRI-1-1109 Crease Nail Puller
For the removal of nails when worn down inside the grooves of the shoe, 12" s/s



BRI-1-1111 Shoe Spreader
The design of this tool is simply a levering device with a differential length of arms on the head. The shoe there by opening it exactly where you want.



BRI-1-1113 Clinch Cutter Powder Coated
Clinch cutter on one end and pritchel on opposite end, Powder coated c/s.



BRI-1-1114 Farrier Hammer

Grooming Equipments



BRI-1-1115 Horse Machine

Removes unwanted whiskers and hairs. American pattern, Sand caste, Chromed.



BRI-1-1116 Horse Machine

Removes unwanted hair and whiskers from Horse & Cattle, Putor, German pattern. Red wooden handle.



BRI-1-1117 Machine For Sheep

Removes unwanted hair from Sheep Putor, German pattern, Red wooden handle.



BRI-1-1118 Hand Clipper

Spring type sheep shear (For Horse, Camel and Donkey)

BRI-1-1119



For sheep



BRI-1-1100 Bot Egg Knife Plastic

The safe efficient way to remove bot eggs. It has no sharp edges to gouge the horse and can be used by the youngest rider with safety. This bot egg knife is 8.75" long with a stainless steel blade.



BRI-1-1100 Bot Egg Knife Wooden

This bot egg knife is 8.75" long with a stainless steel blade and durable wooden handle.



BRI-1-1100 Fetlock Shears

Can be used in either hands with blades curving left or right. Stainless steel.



BRI-1-1100 Post Mortem Scissor

Heavy-duty scissors with one serrated edge. Sand caste, Stainless steel, length 7-3/4". Also available with black finished rings.

Hoof Rasps



BRI-1-1100 Flat Rasp

A flat rasp with rasp cut on one side and double cut file on other. Super sharp. Length 36cm c/s.



BRI-1-1100 Flat Rasp

A flat rasp with rasp cut on one side and double cut file on other. Super sharp. Length 36cm c/s.



BRI-1-1100 Hoof Rasp X Style

Rasp teeth on one side with file teeth on other. Length 31cm c/s.



BRI-1-1100 Hoof Rasp X Style

Rasp teeth on one side with file teeth on other. Length 31cm c/s.



BRI-1-1100 Hoof Rasp

Rasp teeth on two side with file teeth on other. Length 31cm c/s.



BRI-1-1100 4-Way Hoof Rasp

The file has different filing grades: fine, medium fine, medium. This allows you to file the hooves to proper length while shaping and smoothing the rough and sharp edges. The 4-Way Hoof Rasp is the best tool for smoothing all areas of the animal's hooves. Length 35cm c/s.

Hoof & Claw Instruments



BRI-1-1130 Hoof Tester

Test hoof hardness with this popular, steel hoof tester with scored edges for better gripping. Available sizes 13", 13.25", 16" s/s.



BRI-1-1131 Hoof Tester Heavy Duty

Available sizes 13", 13.5" s/s.



BRI-1-1132 Hoof Tester

Length 16" s/s.



BRI-1-1133 Large Hoof Tester

The large model has 12" reins, is 19" long and opens to 9.5". Stainless steel with pointed teeth.



BRI-1-1134 Small Hoof Tester

The small model has 9" reins, is 15" long and opens to 7". Stainless steel with pointed teeth.



BRI-1-1135 Large Modified Hoof Tester

The modified model has 12" reins, is 19" long and opens to 9.5". Stainless steel with reduced teeth



BRI-1-1136 Hoof Trimmer

Durable steel with bolted construction. Green protective coating covers the steel. 1-1/4" opening between trimmer blades with double compound action. Sharp blades are firmly attached with three recessed screws. Black rubber grip handles. Length 17" c/s.



BRI-1-1137 Nail Clipper

A nail clipper with sure-grip handles for cutting dogs nails. Stainless steel blades are designed to give precision cutting action. Flip-open safety lock protects cutting edge when not in use. 6-3/8" long.


Hoof & Claw Instruments



BRI-1-1138 Foot Rot Hoof Shear Carbon Steel
Carbon steel handles: Length 9.5".



BRI-1-1139 Foot Rot Hoof Shear Aluminum
Ideal for both sheep and goat. Length 9.5".



BRI-1-1140 Foot Rot Hoof Shear Serrated
Primarily used for sheep: Length 9" s/s.



BRI-1-1141 Foot Rot Hoof Shear Plane
Available with plane blades and in leather holster 9" s/s.



BRI-1-1142 Trimming Shears
Red enameled handles with extra-large double bow Length 14".



BRI-1-1143 Single Bow Economy Sheep Shears
6" shearing blades. Overall Length 12".



BRI-1-1144 Hoof Pick-Vinyl Coated Handle
Sturdy, vinyl coated handle: Length 5.5" s/s.



BRI-1-1145 Hoof Pick Large - Vinyl Coated Handle
Pick with vinyl coated handle, s/s.



BRI-1-1147 Hoof Knife Kit
Available in compact, elegant leather case with zipper. Set includes Wooden handle with 4 different types of Interchangeable blades.




BRI-1-1146 Hoof Pick
Chrome plated steel rod hoof pick with a screw driver end which can be used to tighten Chicago screws.



BRI-1-1148 Hoof Knife Kit
Available in compact, elegant leather case with zipper. Set includes Wooden handle with 6 different types of Interchangeable blades.

Hoof & Claw Instruments




BRI-1-1149 Dog Nail Scissor

For cutting dog and other animals nails. Length 4.5", s/s.



BRI-1-1150 Cat Nail Scissor

For cutting cat nails. Length 3.5", s/s.



BRI-1-1151 Swiss Hoof Knife

Used for cutting and cleaning large areas of the sole in cattle and other large animals. Overall length: 7" s/s. Polished wooden grip.

BRI-1-1152 Hoof Knife

Double edged. Length 7.5" s/s. Polished wooden grip.



BRI-1-1153 Hoof Knife Small Size

Double edged. Length 6.5" s/s. Polished wooden grip



BRI-1-1154 Oval Hoof Knife

Available in wooden handle, can be used left or right handed. Extra long shaft makes it easier and more accurate to use



BRI-1-1155 Oval Hoof Knife

Available in wooden handle, can be used left or right handed.



Hoof & Claw Instruments



Right

BRI-1-1156

Left

Hoof Knife Regular

Regular blade. This durable blade is held by two rivets in the hardwood handle and it extends 2" into handle.
1/2" blade at heel, 3" long. 8.5" long overall. Wooden grip. Available in left and right handed grip.



Right

BRI-1-1157

Left

Hoof Knife

Double edged. Polished wooden grip.
Available in left and right handed grip.



Right

BRI-1-1158

Left

Hoof Knife

Polished wooden grip. Available
in left and right handed grip.



Right

BRI-1-1159

Left

Hoof Knife

Narrow blade, Polished wooden grip.
Available in left and right handed grip.



Right

BRI-1-1160

Left

Hoof Knife

Double edged. Polished wooden grip.
Available in left and right handed grip.



Right

BRI-1-1161

Left

Hoof Knife

Curved blade. Polished wooden grip.
Available in left and right handed grip



Right

BRI-1-1162

Left

Hoof Knife

Extra long shaft makes it easier and more accurate to use.
Polished wooden grip. Available in left and right handed grip.



Right

BRI-1-1163

Left

Hoof Knife

Narrow blade. Fine wooden grip.
Available in left and right handed grip.

Hoof & Claw Instruments



Right

BRI-1-1164

Left

Hoof Knife

Wide blade. Fine wooden grip.
Available in left and right handed grip.



Right

BRI-1-1165

Left

Hoof Knife

Regular blade. Plastic grip.
Available in left and right handed grip.



Right

BRI-1-1166

Left

Hoof Knife

Double edged. Plastic grip.
Available in left and right handed grip.



Right

BRI-1-1167

Left

Hoof Knife

Double edged. Plastic grip.
Available in left and right handed grip.



Right

BRI-1-1168

Left

Hoof Knife

Narrow blade. Plastic grip.
Available in left and right handed grip.



Right

BRI-1-1169

Left

Hoof Knife

Narrow blade. Double edged. Plastic grip.
Available in left and right handed grip



Right

BRI-1-1170

Left

Hoof Knife Backolyte Grip


Narrow blade. Backolyte grip.
Available in left and right handed grip.

Injectors & Extractors



BRI-1-1171 Europlex Syringe

Europlex metal gun syringe 30ml, 50ml sizes
(Also available with casing).



BRI-1-1172 Europlex Syringe

Europlex metal gun syringe 30ml, 50ml sizes
(Also available with casing).



BRI-1-1175 Automatic Syringe

Automatic Syringe
Available sizes: 5ml, 10ml, 15ml




BRI-1-1173 Bloat Trocar

The flat-to-round shape has made this wooden-handled trocar very popular for venting animal gas accumulations. Includes stainless steel cannula and protector cap. Cannula length is 6". Also available in 5" length.




BRI-1-1174 Cattle Trocar

This cattle trocar is used to relieve bloating. Complete with cannula, a hollow tube which remains when trocar is withdrawn to allow gas to escape. Stainless point. Metal handle 9.5" L. Tube is 5-1/2" L. 3/8" I.D.




BRI-1-1176 Automatic Syringe

Automatic Syringe
Available sizes: 30ml



BRI-1-1177 Automatic Syringe

Automatic Syringe
Available sizes: 60ml



BRI-1-1178 Semi Auto Syringe

Semi Auto Syringe
Available sizes: 70ml



BRI-1-1179 Semi Auto Syringe

Semi Auto Syringe
Available sizes: 150ml

Obstetrical Instruments



BRI-1-1180 Krey Hook
Length 8" s/s.



BRI-1-1181 Linde's Embryotomy Knife
Length 5" s/s



BRI-1-1181 Ostertag's Eye Hook Blunt
Length 2", Hole Diameter 5/8"
and Hook Length 3/4" s/s.



BRI-1-1183 Moore's OB Chain Handle
Stainless steel OB handle with one- piece
construction. 5.5"



BRI-1-1184 Moore's OB Chain Handle
With Plastic Grip Handle aids in use of
chains, 5.5" s/s.



BRI-1-1185 OB Wire Guide
The instrument will bring your wire
around any tight or difficult location
on the fetus. Length 7" s/s.




BRI-1-1186 OB Chain
Heavy duty obstetrical chains assist
in the delivery of the calves.
Length 60" s/s.




BRI-1-1187 OB Chain
Heavy duty obstetrical chains assist
in the delivery of the calves.
Length 30" s/s.



BRI-1-1188 Swine Cable OB Snare
This durable, swine cable obstetrical snare
is made from aircraft cable.
Overall length 18" s/s.



BRI-1-1189 OB Wire Guide
The instrument will bring your wire
around any tight or difficult location
on the fetus. Length 9" s/s.



BRI-1-1191 Farrowing Forceps
Farrowing Forceps For Piglets. Length 20.5" s/s.



BRI-1-1190 Sutter OB Handle
Knurled non-slip grip. Length 4" c/s.



BRI-1-1192 Farrowing Snare
Shaped end allows operator to direct head after cable
is placed behind ears or legs. Length 15" c/s.


Plaster Instruments



BRI-1-1193 Cast Spreader
Cast spreaders with 5cm wide blades.
Available sizes 10.5", 11" s/s.



BRI-1-1194 Cast Spreaders (Henning)
3-Pong cast spreader 9". Spring-loaded handle for
one-hand operation. cutting beyond cast thickness.




BRI-1-1195 Plaster Shears
Available size 17.25" s/s.



BRI-1-1196 Plaster Shears
Plaster shears with serrated blades.



BRI-1-1197 Plaster Shears (Brun/Bohler)
They can cut through cloth, plastic and metal. They have serrated blades for non-slip cutting. Stainless steel. Length 10"



BRI-1-1198 Plaster Shears (Stille)
Plaster shears with serrated blades.
Available sizes 9", 10" 15" s/s.



BRI-1-1199 Ring Cutter 16 cm



BRI-1-1200 Universal Scissors
It can cut through almost anything: a multitude of uses around any clinic. Autoclavable.
Stainless steel: Length 7.5", 5.5".

Restraining Instruments



BRI-1-1203 Premium Hog Holder With Lock

An improved hog holder made of 3/16" galvanized aircraft cable with a cable stop that avoids kicking. Durable, high-strength tubing of carbon steel, nickel-plated.



BRI-1-1204 Ring-O-Matic Hog Catcher

Snares and holds any weight hog. Loop is constructed of galvanized aircraft cable and the handle is coated with a rust-resisting finish. Handle is 22" in length with a 20" loop of aircraft cable.



BRI-1-1205 Shorty Hog Holder

Easy, quick way to control pigs in tight spaces. Comfort-grip rubber handle with galvanized aircraft cable loop. Approximately 8" long.



BRI-1-1206 Canine Mouth Gag

Spring action with nylon inserts for less trauma. Three sizes available: Small: 4.5" s/s, Medium: 5.75" s/s, Large: 6.25" s/s.

R-1615



BRI-1-1207 Cattle Twitch / Nose Twitch

For help in restraining horses or cattle. Nylon coated handles for firm grip. Length 17" s/s.



BRI-1-1208 Bull Holder

Stainless steel sand casted. PVC grip. Available sizes 14", 15" s/s.



BRI-1-1209 Bull Lead With Chain

Place stainless steel lead into animal's nostrils and close gently. Secure handles with heavy-duty 13 1/2" chain. Also available without chain.



BRI-1-1210 Harm's Bull Holder

Made of Stainless steel.



BRI-1-1211 Cattle Leader Aluminum

The ideal type lead for show animals - gives the nose a shorter appearance than other type leads. Made of strong, Permits animal to eat without removal. 2-1/2" x 1-1/4" x 3-3/4" L.



BRI-1-1212 Cattle Leader Brass

The ideal type lead for show animals - gives the nose a shorter appearance than other type leads. Permits animal to eat without removal.



BRI-1-1213 Heavy Duty Bull Lead


With coiled spring and protective release sleeve. Made of Stainless Steel.



BRI-1-1214 Kantslip Cattle Lead

Made with fine quality holding spring. Opens to approximately 2-5/8". Overall width 3", overall length 5" brass, c/p.

Ring Applicators



BRI-1-1225 Snout Pig Ring Applicator

Sand cast, Stainless steel snout pig ring applicator. Length 6" s/s.




BRI-1-1227 Bull Nose Punch

For piercing the bull's nose cartilage prior to insertion of the nose ring. This tool makes the job more animal friendly.



BRI-1-1229 Bull Ring Applicator Heavy Duty

Stainless Steel bull ring applicator. Excellent Quality. This is a far superior product than the standard steel applicators.



BRI-1-1226 Nose Ring Applicator

Make pig ringing quick, easy and painless. Rings go between nostrils. Length 10".



BRI-1-1228 Bull Ring Applicator Economy

This metal applicator will do rings up to 77mm x 10mm (3 inch) Length overall is 11".



BRI-1-1230 Bull Ring Stainless

The special sharp cutting points pierce the nasal septum easily with no pre-incision. Size's Available: 2", 2.5", 3", 3.5".



BRI-1-1231 Bull Ring Brass

Brass bull ring is self-piercing and extremely durable. A long, sharp piercing point ensures a clean cut. Size's Available: 2", 2.5" 3", 3.5"

Teat Instruments



BRI-1-1232 Teat Slitter

Opens test sphincter for greater milk flow. Available with loops. Length 8", brass and s/s.



BRI-1-1234 Hug's Teat Tumor Extractor

Effective tool used to remove fibrous tissue to prevent calculi from forming: Length 6", brass and s/s.



BRI-1-1236 Cornell teat curette

Stainless steel with spiral cutting edge 6"




BRI-1-1238 Hug's Teat Knife / Teat Lancet

Length 6", blade 0.5", s/s.



BRI-1-1239 Sudan Blood Taking Needle

Length 3.5", brass and s/s.



BRI-1-1233 Teat Tumor Extractor

Effective tool used to remove fibrous tissue to prevent calculi from forming. Length 7.5", brass and s/s.



BRI-1-1235 Danish Teat Slitter

Extra thin model. Blades are opened by loosening the small burr on side of the instrument: Length 7.5", brass and s/s.



BRI-1-1237 Teat Syphons


Milking Tube Length 3.75", brass and s/s.



BRI-1-1240 Milking Tube


Length 4.25", brass and s/s.

Vaginal Specula



BRI-1-1224 Polansky

3 Blades for wide exposure. Butterfly screw with 2 blades.
Blades Length 27cm, Overall Length 35 cm, Max.
opening 10 cm, s/s.



BRI-1-1224 Polansky Junior

3 blades for wide exposure. Butterfly
with 2 blades Blades Length 24cm, s/s.



BRI-1-1224 Winterton

Blades Length 10cm, s/s.



BRI-1-1224 Russian

Blades Length 23cm, s/s.



BRI-1-1224 Graves

Blades Length 15cm, s/s.

Snaffles. Dee. Ring Snaffle Bits



BRI-2-1246 Rings stainless steel

4" to 5-1/2" mouth size available. Ring varies from 45mm to 75mm.



BRI-2-1247 Rings stainless steel

Mouth and rings made of stainless steel.



BRI-2-1248 Rings stainless steel

Rings stainless steel. Mouth german steel.



BRI-2-1249 Rings stainless steel

Rings stainless steel. Mouth carbon steel with leather cover.



BRI-2-1250 Rings stainless steel

Stainless steel rings. Mouth made of rubber.



BRI-2-1251 Rings stainless steel

D bit with 5" rubber covered mouth. Stainless steel. Also available in chrome plating.




BRI-2-1252 D ring snaffle

made of stainless steel. 4-1/2" to 5-1/2" mouth sizes available.



BRI-2-1253 Rings stainless steel

Stainless steel. mouth made of german steel including chain.



BRI-2-1254 Rings stainless steel

Rings and mouth made of stainless steel.



BRI-2-1255 Ring Stainless steel

Stainless steel rings. Link and mouth made of stainless steel.



BRI-2-1256 Ring Stainless steel

Stainless steel rings. Link made of red copper. Stainless steel mouth.



BRI-2-1257 Ring Stainless steel

Stainless steel rings. Link and mouth made of red copper.



BRI-2-1258 Ring Stainless steel

Brass rings. Link and mouth made of brass.



BRI-2-1259 Ring Stainless steel

Stainless steel.



BRI-2-1260 Ring Stainless steel

Stainless steel. Mouth made of german steel.



BRI-2-1261 Ring Stainless steel

Stainless steel.



BRI-2-1262 Ring Stainless steel

Stainless steel.



BRI-2-1263 Ring Stainless steel

Stainless steel. Mouth german steel including chain.



BRI-2-1264 Ring Stainless steel

Stainless steel.



BRI-2-1265 Ring Stainless steel

Link german silver. Mouth made of stainless steel.



BRI-2-1266 Ring Stainless steel

Stainless steel pelham including chain.


Snaffles. Dee. Ring Snaffle Bits




BRI-2-1267 Snaffle Bits
Stainless steel including chain.




BRI-2-1268 Snaffle Bits
Stainless steel. Link German steel including chain.



BRI-2-1269 Snaffle Bits
Rings stainless steel. Mouth German steel.



BRI-2-1270 Snaffle Bits
Rings and mouth made of stainless steel.



BRI-2-1271 Snaffle Bits
Stainless steel including chain.



BRI-2-1272 Snaffle Bits
Rings stainless steel. Rubber covered mouth.



BRI-2-1273 Snaffle Bits
Stainless steel including chain.



BRI-2-1274 Snaffle Bits
Stainless steel hollow mouth including chain.




BRI-2-1275 Snaffle Bits
Stainless steel.




BRI-2-1276 Snaffle Bits
Stainless steel.



BRI-2-1277 Snaffle Bits
Full cheek with stainless steel mouth.



BRI-2-1278 Snaffle Bits
Stainless steel.



BRI-2-1279 Snaffle Bits
Stainless steel.

Spurs and Stirrups



BRI-2-1280 Spurs Stainless Steel

Stainless steel spurs with brass buttons and rowel. Also available in chrome plating. 2pcs set.



BRI-2-1281 Spurs Stainless Steel

Stainless steel spurs. 2pcs set.



BRI-2-1282 Spurs Stainless Steel

Classic hunt and jump spurs in the popular loop end style. Stainless steel. 2pcs set.



BRI-2-1283 Spurs Stainless Steel

Fine engraving with loop end style. Stainless steel. 2pcs set.



BRI-2-1284 Spurs Stainless Steel

Stainless steel spurs. Gold plated. 2pcs set.



BRI-2-1285 Spurs Stainless Steel

Carbon steel spurs with powder coating. Buttons and rowel made of brass.



BRI-2-1286 Spurs Stainless Steel

Carbon steel spurs powder coated with brass buttons. 16 point rowel. 2pcs set.



BRI-2-1287 Spurs Stainless Steel

Brass spurs with 10 point rowel. 2pcs set.



BRI-2-1288 German steel stirrup

Including rubber, 2pcs set. Sizes Available: 4", 4 1/4", 4 1/2", 4 3/4", 5".



BRI-2-1289 Stainless steel stirrup

Including rubber, 2pcs set. Sizes Available: 4", 4 1/4", 4 1/2", 4 3/4", 5".



BRI-2-1290 Stirrup

Stainless steel stirrup, flexible. Including rubber, 2pcs set.



BRI-2-1291 Stainless steel stirrup

Including rubber, 2pcs set. Sizes Available: 4", 4 1/4", 4 1/2", 4 3/4", 5".



BRI-2-1292 Stainless steel stirrup

Including rubber, 2pcs set. Sizes Available: 4", 4 1/4", 4 1/2", 4 3/4", 5".



BRI-2-1293 Stainless steel stirrup

Stainless steel stirrup, Including rubber, 2pcs set. Sizes Available: 4", 4 1/4", 4 1/2", 4 3/4", 5".



BRI-2-1294 Stainless steel stirrup

Stainless steel stirrup, Including rubber, 2pcs set. Sizes Available: 4", 4 1/4", 4 1/2", 4 3/4", 5".



BRI-2-1295 Stainless steel stirrup







Especially designed for polo players. Sizes Available: 4", 4 1/4", 4 1/2", 4 3/4", 5".



BRI-2-1296 Rubber for stirrups

Rubber for stirrups. 2pcs set.

Surgical Forceps



Towel Clamps



Dressing and Tissue Forceps

BRI-3-1392

- (a) Adson Serrated Forceps **Size:** 12.5cm
- (b) Adson Brown Dissecting Forceps **Size:** 12.5cm
- (c) Adson Dissecting Forceps 1/2 Teeth **Size:** 12.5cm



BRI-3-1393

Lanes Dissecting Forceps 1/2 Teeth
Size: 15cm



BRI-3-1394

Surgical Dressing Forceps
Sizes: 13cm, 14cm, 16cm, 18cm, 20cm, 25cm, 30cm



BRI-3-1395

Standard Tissue Forceps 1/2 Teeth
Sizes: 13cm, 14cm, 16cm, 18cm, 20cm, 25cm, 30cm



BRI-3-1396

Russian Tissue Forceps
Sizes: 16cm, 20cm



Hollowares



BRI-3-1397

Kidney Dishes
Sizes: 6", 8", 10" 12"



BRI-3-1398

Instrument Steel Boxes Box with Lid
250mm x 120mm x 50mm
300mm x 150mm x 50mm
400mm x 200mm x 70mm



BRI-3-1399

Lotion Bowls
200mm x 100mm x 40mm
Sizes: 6", 8", 10"

Scalpel Handles

BRI-3-1400

Scalpel Handle No. 3








BRI-3-1401




Scalpel Handle No. 4



Needle Holders




Scissors





BRI-3-1416 Balfour Abdominal Retractor

Balfour Abdominal Retractor 180mm Opening
Balfour Baby Abdominal Retractor 90mm Opening



BRI-3-1419 Finochietto Rib Retractor

Finochietto Rib Retractor Small 75mm Spread
Finochietto Rib Retractor Large 140mm Spread



BRI-3-1421 Rib spreader Stainless steel

Rib spreader Stainless steel Standard Size,
blades 48mm deep x 65mm wide, 12" Spread, Curved arms



BRI-3-1423 Morris Single Ended Retractors

Morris Single Ended Retractor 5.2x3.8cm
Morris Single Ended Retractor 5.2x5.2cm
Morris Single Ended Retractor 5.2x6.4cm



BRI-3-1417 Gelpi Retractors

Mini Gelpi Retractor 90cm Gelpi Retractor 14cm 5.5"
Gelpi Retractor 18cm 7" Weitlaner Retractor 6.5" Sharp Teet



BRI-3-1418 Travers Retractor

Travers 4x5 Teeth Retractor 20cm 8"



BRI-3-1420 Weitlaner Retractors

Weitlaner Retractor 7" Blunt Teeth Small Weitlaner
Retractor 5.5" Sharp 3x4 Teeth Weitlaner Retractor 6.5" Blunt 3x4 Teeth



BRI-3-1422 Senn Miller Retractors

Senn Miller Retractor Blunt Teeth
Senn Miller Retractor Sharp Teeth



BRI-3-1424 Langenbeck Retractors

Langenbeck Retractors 2.3x0.7cm
Langenbeck Retractors 4.4x1.3cm
Langenbeck Retractors 6.5x2.5cm
Langenbeck Retractors 9.0x2.0cm

Equine Surgical Instruments




BRI-3-1425

Cervix Forceps Equine Realted Size 21" 53cm Easy to hold pistol handle with step locking ratchet mechanism and easy to release one click hook for ease of operation.




BRI-3-1426

Alligator Forceps 14cm 5.5" Serrated



BRI-3-1427

Jackson uterine biopsy forceps
Stainless steel biopsy forceps has a 4mm x 28mm cutting area and 60cm long shaft. Popular for equine uterine biopsy



BRI-3-1428

Mare Biopsy Punch 26". Mare Uterine Biopsy Punch is an equine reproduction product for ues in artificial insemination in horses. Mare biopsy punch with pistol grip handle and locking nut. (26") Long

Miscellaneous Instruments



BRI-3-1429

Equine Stomach Pump
Size: 15.5 (Small) and 21(Large).
Supplied with 10 feet long tube



BRI-3-1430

Mini fibre optic Otoscope with Pouch



BRI-3-1431

Mole Trap



BRI-3-1432

Digital Thermometer C or F
(Picture for illustration Purposes only)



BRI-3-1433

Embryotome Thygesen 85 cm



BRI-3-1434

Ear Syringe Complete:
50ml, 75ml, 100ml, 150ml, 200ml

Root Elevators



BRI-3-1435

Winged Elevators Standard Stubby 1mm-6mm



BRI-3-1436

Winged Elevators Standard Handle 1mm-6mm



BRI-3-1437

Straight Root Tip Elevator Standard Handle



BRI-3-1438

Periosteal Elevator Standard Handle



BRI-3-1439

Straight Root Tip Stubby Handle Elevators 1mm-2mm(Also Set of 2)



BRI-3-1440

Set of Winged Tip Elevators in sizes 1-4mm




BRI-3-1441

Set of Winged Tip Elevators in sizes 1-4mm

Root Elevators

BRI-3-1442
Canine Root Elevator



BRI-3-1443
Cryer Elevator Set of 2



BRI-3-1444
Coupland Dental Elevator Set of 1mm-3mm(Set of 3)



BRI-3-1445
Elevator Sterilization Tray



BRI-3-1446
Elevator Sterilization Tray

Scalers & Probes for Cats & Dogs



BRI-3-1447
Explorer/Measuring Probe Standard



BRI-3-1449
Jaquette Supragingival Scaler Double Ended



BRI-3-1448
Periosteal Elevator Double Ended



BRI-3-1450
Subgingival Curette



BRI-3-1451
Sickle Scaler Normal

Specialized Dog Scalers



BRI-3-1452
McCall Dog Scaler, Dual Ended



BRI-3-1454
Dog Scaler Left



BRI-3-1456
Periosteal Elevator No 9



BRI-3-1453
Kaplan Dog Scaler



BRI-3-1455
Dog Scaler Right



BRI-3-1457
Examination Mirror Single Ended used for Small Animal Dentistry

Dental Extraction Forceps & Mouth Retractors



BRI-3-1458
Extraction Forceps Angled Jaw, 12 cm



BRI-3-1459
Tooth Extraction Forceps, Angled & Straight Jaw, 2.5 cm



BRI-3-1460
Small Animal Forceps Fragment, 16.5 cm



BRI-3-1461
Root Tip Extraction Forcep, 15 cm



BRI-3-1462
Cheek Retractor Minnesota



BRI-3-1463
Small Animal Mouth Gag Openers (Set of 5 Gags)

- Mouth Gag for Dog 16 cm
- Mouth Gag for Dog 12.7 cm
- Cat/Dog 8.8 cm
- Cat/Dog 5.5 cm
- Special Cat 2.7 cm

VETERINARY INSTRUMENTS CATALOGUE

Ask For it!

info@bergmed.com

*Our Products Are Guaranteed And Meet The Relevant
Worldwide Standards And Our Products Are
Manufactured Entirely Under In Factory.*



+92 (308) 160 7246



www.bergmed.com



bergmedint@gmail.com | info@bergmed.com



Small industrial Area, Sialkot - 51310 Pakistan